NOTICE OF APPROVAL

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Issued by statutory authority of the Minister of Industry
for:

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Chromatographe en phase gazeuse

APPLICANT

Norwest Soil Research Ltd.
9958 - 67 Avenue
Edmonton, Alberta
T6E 0P5

REQUÉRANT

MANUFACTURER

Agilent Technologies, Inc.
2850 Centerville Road
Wilmington, Delaware
19808-1610

FABRICANT

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Agilent 6890A (G1540A) Gas
Chromatograph/
Chromatographe en phase gazeuse
Agilent 6890A (G1540A)
Hewlett Packard 5890

RATING/ CLASSEMENT

20°C to/à 27 °C

Helium and Hydrogen only/seulement

NOTE: This approval applies only to meters, the design, composition, construction and performance of which are, in every material respect, identical to that described in the material submitted, and that are typified by samples submitted by the applicant for evaluation for approval in accordance with sections 13 and 14 of the Electricity and Gas Inspection Regulations. The following is a summary of the principal features only.

SUMMARY DESCRIPTION:

The Agilent 6890A (G1540A) Series Gas Chromatograph analyzes the composition of natural gas. The gas chromatograph separates the following components of natural gas:

S	oxygen
S	nitrogen
S	methane
S	carbon dioxide
S	ethane
S	propane
S	iso-butane
S	n-butane
S	iso-pentane
S	n-pentane
S	hexanes
S	heptanes
S	octanes
S	nonanes
S	decanes

The individual component concentrations are used to calculate the relative density and heating value of the sample. The detection of components is performed by a TCD (thermal conductivity detector) and a FID (flame ionization detector) which provide an output signal proportional to the concentration.

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Natural gas is injected into the sampling valve and is carried by the carrier gas, helium through four columns, inside the oven. The different components move through the columns at different rates, thereby allowing identification based on the retention time of each component.

The operator interface consists of a display, status board, and keyboard on the face of the gas chromatograph as well as a user interface program (Chemstation software) running on a personal computer. The Chemstation software analyses the data and displays the chromatogram. The Chemstation software version number can be viewed by selecting the Help menu and then "About" in the Chemstation software program.

A Foxpro software application entitled Analysis is used to determine calorific power and relative density, using GPA 2172 and GPA 2145 at 15°C and 101.325 kPa. The Foxpro Analysis program version is identified by date. The date can be viewed by opening the Windows Explorer application and looking beside the Foxpro Analysis program file entitled "anal_sys.exe" located in the "c:\an_prog" directory.

The Agilent 6890A gas chromatograph's firmware version can be viewed by powering down and then powering up again the gas chromatograph. The firmware version will appear on the screen during boot up.

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A Hewlett Packard model 5890 analyzer can be used to improve the accuracy of the gas analysis and hence the relative density and heating value calculations. Because the 6890 A uses helium as a carrier gas it can not be used to detect helium so detecting helium (and hydrogen) with the 5890 and then feeding that result into the Foxpro software along with the output from the model 6890 A enhances the analysis by making the heating value and the relative density figures more reliable. The model 5890 is used only to measure these two gases and only when they are anticipated components in the gas sample under analysis.

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Ambient
Temp. Range: 20°C to 27°C

Oven Temp. Range: - 80°C to 450°C

Column 1: 6' x 1/8" - 35% DC200/500
on Chromosorb P/AW
- 80/100 Back Flush Packed

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Foxpro, Analysis Software,
March 29, 2001, June 10th
2002

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Gas Chromatograph

Adjunct: HP 5890
Model No. 5890A (19234B)
Serial No. 2413AO 2428
NWL Inv. No. 1985
Voltage: 120VAC, +5% / -10%
Frequency: 47.5 to 67 Hz
Power consumption: 2200VA max
Firmware/software Div 43 Rev D 5890-6000
Carrier Gas: argon
Columns: mole seive 13 x 9' 80% mesh

Detector: thermal conductivity (TCD)

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MARKING REQUIREMENTS

Marking shall be in accordance with Section 3.5 and
20-3.1 of LMB-EG-08.

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REVISIONS**AG-0489 Rev. 1**

Revision 1 approves the Hewlett Packard model 5890 chromatograph, calibrated for helium and hydrogen only, to enhance compositional resolution of the model 6890 A chromatograph.

EVALUATED BY**AG-0489**

Judy Farwick
Complex Approvals Examiner
Tel: (613) 946-8185
Fax: (613) 952-1754

AG-0489 Revision 1

Graham Collins
Approval Examiner
Tel (613) 941-0605
Fax (613) 952-1754

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
Hewlett Packard 5890

APPROVAL:

The design, composition, construction and performance of the meter type(s) identified herein have been evaluated in accordance with regulations and specifications established under the Electricity and Gas Inspection Act. Approval is hereby granted accordingly pursuant to subsection 9(4) of the said Act.

The sealing, marking, installation, use and manner of use of meters are subject to inspection in accordance with regulations and specifications established under the Electricity and Gas Inspection Act. The sealing and marking requirements are set forth in specifications established pursuant to section 18 of the Electricity and Gas Inspection Regulations. Installation and use requirements are set forth in specifications established pursuant to section 12 of the Regulations. Verification of conformity is required in addition to this approval for all metering devices excepting instrument transformers. Inquiries regarding inspection and verification should be addressed to the local inspection office of Industry Canada.

Original copy signed by:

René Magnan, P. Eng
Director
Approval Services Laboratory

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RÉVISIONS**AG-0489 Révision 1**
ÉVALUÉ PAR**AG-0489**
AG-0489 Révision 1