**NOTICE OF APPROVAL**

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Issued by statutory authority of the Minister of Industry  
for:

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**REQUÉRANT**

**MANUFACTURER**

Varian Chromatography Systems  
2700 Mitchell Drive  
Walnut Creek California  
USA 94598-1675

**FABRICANT**

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CP-3800 Analyzer "B"

**RATING/ CLASSEMENT**

**NOTE:** This approval applies only to meters, the design, composition, construction and performance of which are, in every material respect, identical to that described in the material submitted, and that are typified by samples submitted by the applicant for evaluation for approval in accordance with sections 13 and 14 of the Electricity and Gas Inspection Regulations. The following is a summary of the principal features only.

### SUMMARY DESCRIPTION:

The Varian CP-3800 Analyzer "B" gas chromatograph analyses the gas sample and reports the gas composition. This information is then entered into the AltaGas Utilities Data 01 excel spreadsheet which calculates the heating value and relative density of the gas sample.

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#### Valves

- Stream selection
  - Valco model EMTCS10UW
  - 10-port high pressure valve
- Gas sampling
  - Valco model EMTCS10UW
  - 10-port high pressure valve
- Back-Flush
  - Valco model EMT2CSD12MWE
  - 12-port valve

#### Operating Gas

- Carrier: Helium, 99.999% purity.
- Make-up: Helium, 99.999% purity
- FID detector fuel gas
  - Hydrogen, 99.999% purity
  - Ultra Zero Grade Air

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## Metrological Functions

### Gas Composition

The Varian CP 3800 is approved for determining relative concentrations of the following components in a gas sample:

- N<sub>2</sub>, nitrogen
- O<sub>2</sub>, oxygen
- CO<sub>2</sub>, carbon dioxide
- CH<sub>4</sub>, methane
- C<sub>2</sub>H<sub>6</sub>, ethane
- C<sub>3</sub>H<sub>8</sub>, propane
- i-C<sub>4</sub>H<sub>10</sub>, iso-butane
- n-C<sub>4</sub>H<sub>10</sub>, normal butane
- i-C<sub>5</sub>H<sub>12</sub>, iso-pentane
- n-C<sub>5</sub>H<sub>12</sub>, normal pentane
- C<sub>6</sub>H<sub>14</sub>+, hexane and above

### Energy Calculations

The AltaGas Utilities Data 01 Excel Spreadsheet uses the method outlined in AGA-5 for determining the heating value and relative density from the sample gas composition determined by the chromatograph.

## Firmware

The following firmware version is approved for use on the Varian CP 3800:

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- 3.3.5D

**Software**Galaxie Workstation

The following version of the “Galaxie Workstation” are approved for use:

- 1.8.508.1

A PC running “Galaxie Workstation” software can be used to communicate with the model 3800 gas chromatograph. Use of this software permits configuration, calibration and data retrieval.

AltaGas Utilities Data 01 Excel Spreadsheet

The following version of the Data 01 spreadsheet is approved:

- Rev. 95-03-07

**SPECIFICATIONS**

- Operating temperature range +10° to +40°C
- Power supply input voltage 120 Vac, 60 Hz
- Power supply input power 2400 VA

**MARKINGS**

Marking shall be in accordance with LMB-EG-08 sections:

- 3-5.1 and
- 20-3.1.

**SEALING PROVISIONS**

Sealing provisions are not required for gas chromatographs.

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**VERIFICATION**Energy Calculations

The true calorific power and relative density of the sample gas shall be calculated using the true relative concentrations of all components of the sample gas as established in accordance with the Gas Processors Association standards, *Calculation of Gross Heating Value, Relative Density and Compressibility Factor for Natural Gas Mixtures from Compositional Analysis*, GPA 2172, and *Table of Physical Constants for Hydrocarbons and Other Compounds of Interest to the Natural Industry*, GPA 2145.

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**FIGURES**

**APPROVAL:**

The design, composition, construction and performance of the meter type(s) identified herein have been evaluated in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. Approval is hereby granted accordingly pursuant to subsection 9(4) of the said Act.

The sealing, marking, installation, use and manner of use of meters are subject to inspection in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. The sealing and marking requirements are set forth in specifications established pursuant to section 18 of the Electricity and Gas Inspection Regulations. Installation and use requirements are set forth in specifications established pursuant to section 12 of the Regulations. Verification of conformity is required in addition to this approval for all metering devices excepting instrument transformers. Inquiries regarding inspection and verification should be addressed to the local office of Measurement Canada.

Original signed by:

Patrick J. Hardock, P.Eng.  
Senior Engineer – Gas Measurement  
Engineering and Laboratory Services Directorate

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