**NOTICE OF APPROVAL**

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Issued by statutory authority of the <sup>Mesures Canada</sup> Minister of Industry for:

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Chromatographe en phase gazeuse

**APPLICANT**

ORTECH Environment  
804 Southdown Road  
Mississauga, On  
Canada, L5J 2Y4

**REQUÉRANT**

**MANUFACTURER**

Hewlet Packard  
3000 Haover Street  
Palo Alta, California  
USA, 94304

**FABRICANT**

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5890 Series II - 10643

**RATING/ CLASSEMENT**

Heating Value/Pouvoir calorifique : de 400 to/à 1 200 BTU/ft<sup>3</sup>/pi<sup>3</sup>  
Relative Density/Densité relative : de 0,500 to/à 1,035

**NOTE:** This approval applies only to meters, the design, composition, construction and performance of which are, in every material respect, identical to that described in the material submitted, and that are typified by samples submitted by the applicant for evaluation for approval in accordance with sections 13 and 14 of the Electricity and Gas Inspection Regulations. The following is a summary of the principal features only.

### SUMMARY DESCRIPTION:

The chromatograph consists of two separation columns, one TCD detector, an integrated controller assembly and a PC to determine the mole percent composition of natural gas and landfill gas having a range of energy density from 400 to 1200 Btu/ft<sup>3</sup> and a relative density from 0.500 to 1.035.

The HP 5890 Series II Gas Chromatograph analyses the composition of natural gas and landfill gas.

For natural gas, the sample is analysed using the Chromosorb column. The gas chromatograph separates the following components:

- nitrogen
- oxygen
- methane
- carbon dioxide
- ethane
- propane
- iso-butane
- n-butane
- iso-pentane
- n-pentane
- hexanes+

The results are then compiled using a custom Excel workbook.

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For landfill gas, the sample is analysed using two columns. On the Chromosorb column the following components are separated:

- nitrogen
- methane
- carbon dioxide
- ethane
- propane
- iso-butane
- n-butane
- iso-pentane
- n-pentane
- hexanes+

On the Carbosieve column the following components are separated:

- nitrogen
- carbon dioxide
- oxygen

The combine results from the two runs are the normalize using the raw values to 100% for reporting in mole %.

The results are then compiled using a custom Excel workbook.

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**Carrier Gas:** Helium

### Environment:

- ambient
- 0 - 95% RH (non-condensing)
- indoor / outdoor

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Column #2: 12' x 1/8" Carbosieve SII,  
80/100

Carrier Gas: Helium at 20 to 30 mL/min

Software: HP 3365 Series II Chemstation  
(MII # 11163)

**MARKING REQUIREMENTS**

Markings shall be in accordance with Section 3.5 and 20-3.1 of LMB-EG-08 and section 20-3.1 of S-G-03 (2008-10-08).

Displayed is:

Base temperature and pressure for the calorific power.  
(Also available on printed report).

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
Figure 2 - column / colonne

**APPROVAL:**

The design, composition, construction and performance of the meter type(s) identified herein have been evaluated in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. Approval is hereby granted accordingly pursuant to subsection 9(4) of the said Act.

The sealing, marking, installation, use and manner of use of meters are subject to inspection in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. The sealing and marking requirements are set forth in specifications established pursuant to section 18 of the Electricity and Gas Inspection Regulations. Installation and use requirements are set forth in specifications established pursuant to section 12 of the Regulations. Verification of conformity is required in addition to this approval for all metering devices excepting instrument transformers. Inquiries regarding inspection and verification should be addressed to the local office of Measurement Canada.

**ORIGINAL COPY SIGNED BY:**

Patrick J. Hardock, P.Eng.  
Senior Engineer – Gas Measurement  
Engineering and Laboratory Services Directorate

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