Consumer and
Corporate Affairs Canada

Legal Metrology

Consommation
et Corporations Canada

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APPROVAL No. · N° D'APPROBATION

S.WA-3066 Rev. 2

AUG
AOUT - 3 1988

NOTICE OF APPROVAL

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Issued by statutory authority of the
Minister of Consumer and Corporate
Affairs under application by:

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Fisher Scientific, Division of Allied Canada Inc.
112 Colonnade Road
Nepean, Ontario
K2E 7L6

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NOTE: This approval applies only to
devices, the design, composition,
construction and performance of which
are, in every material respect,
identical to that described in the
information submitted; and are typified
by the sample(s) submitted by the
applicant for evaluation for approval
in accordance with sections 14 and 15
of the Weights and Measures Regula-
tions. The following is a summary of
principal features only.

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Canada

SUMMARY DESCRIPTION:

These are top loading electronic balances. Weight is sensed using an electromagnetic force compensation principle. All models are fitted with levelling feet and a bullseye level. Some models must be operated with a draft shield as the weight increment and sensitivity are such that air currents will prevent the balance from having a true zero. These balances are fitted with + and - indication to the left of the weight display; when both are off, a zero display indicates that the zero position is accurate to $\pm 1/4$ of an increment.

There is only one operator control for these balances which controls, ON/OFF, zero and tare.

- To switch ON - press bar.
- To switch OFF - press bar.
- To zero - press bar. It operates to full scale capacity.
- To tare - press bar. It operates to full scale capacity.

Blanking for over capacity is indicated by a series of flat bars on the weight display. Sealing of the housing is by lead and wire seal through drilled head screws or tamper proof paper seal.

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SUMMARY DESCRIPTION: Continued

All PC models have a Delta Display whereby the last decimal place is blanked out on initial application of weight but as the weight settles it re-appears. This function can be selected by a toggle switch located on the back of the housing. This function speeds up the weighing cycle i.e.; on the PC1616, time for normal weight display is 0.75 sec. to 1.5 sec., with Delta Display ON 0.2 sec. to 0.4 sec.

Note: Delta Range - this is a switch selective mode that increases the increment by a factor of 10 and reduces capacity by the same factor.

These devices are for exclusive use in weighing precious metals or commodities of comparable value and for determining the net content of packages that are not weighed in the presence of the purchaser at the time of sale.

The following accessories are approved and shall be used with the models indicated:

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(1) As these balances have one display only, they must be labelled with instructions that state "THIS BALANCE MUST BE LOCATED SO THAT THE WEIGHT DISPLAY IS VISIBLE TO ALL PARTIES OF ANY TRADE TRANSACTION" (or words to that effect).

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SUMMARY DESCRIPTION: Continued

DESCRIPTION SOMMAIRE: Suite

(2) These balances must be installed on a solid, stable surface free from vibration and in an environment free from air currents.

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(3) These balances are exempted from the following requirements:

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- (a) Zero not to exceed 4% of capacity
- (b) To identify "TARE" when a device is equipped to display only one weight value at a time
- (c) To have levelling feet lockable
- (d) Sections 8, 9.2 and 9.6 of SGM3.

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APPROVAL:

APPROBATION:

The design, composition, construction and performance of the device type(s) identified herein have been evaluated in accordance with regulations and specifications established under the Weights and Measures Act. Approval is granted accordingly, pursuant to subsection 3(1) of the said Act.

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The marking, installation, use and manner of use in trade of devices are subject to inspection in accordance with regulations and specifications established under the Weights and Measures Act. Verification of conformity is required in addition to this approval. Requirements relating to marking are set forth in sections 18 to 26 of the Weights and Measures Regulations. Requirements relating to installation, use and manner of use are set forth in Part V and in specifications established pursuant to section 27 of the said Regulations. Inquiries regarding inspection and certification of conformity should be addressed to the local inspection office of Consumer and Corporate Affairs Canada.

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W.R. Virtue

AUG - 3 1988
AOUT

Chief
Legal Metrology Laboratories

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