Consumer and
Corporate Affairs Canada

Legal Metrology

Consommation
et Corporations Canada

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APPROVAL No. N° D'APPROBATION

S.WA-3331

APR 30 1987

NOTICE OF APPROVAL

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Issued by statutory authority of the
Minister of Consumer and Corporate
Affairs under application by:

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Fisher Scientific
Division of Allied Canada Inc.
112 Colonnade Road
Nepean, Ontario
K2E 7L6

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NOTE: This approval applies only to
devices, the design, composition,
construction and performance of which
are, in every material respect,
identical to that described in the
information submitted; and are typified
by the sample(s) submitted by the
applicant for evaluation for approval
in accordance with sections 14 and 15
of the Weights and Measures Regula-
tions. The following is a summary of
principal features only.

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Canada

SUMMARY DESCRIPTION:

The approved devices are top loading electronic balances. Weight is sensed using an electromagnetic force compensation principle. All models are fitted with lockable levelling feet and a bullseye level.

The devices are equipped with a single operator control which performs the following functions:

To switch ON - press bar.

To switch OFF - lift bar.

To Zero - press bar. It operates to full scale capacity.

To Tare - press bar. It operates to full scale capacity.

A graphic display "Delta Trac" on the right side of the weight display gives an idea of the used and unused portions of the weighing range. When the scale is used, the Delta Trac becomes a plus(+) or minus (-) indicator showing how much weight is on the scale in relation to a target weight.

Blanking for over capacity is indicated by a series of flat bars on the weight display.

Delta Range - this is a mode that decreases the increment size by a factor of 10 and reduces capacity.

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SUMMARY DESCRIPTION: Continued

The models PM and CM incorporate a Delta Display whereby the last decimal place is blanked out on initial application of weight but as the weight settles it re-appears. This function speeds up the weighing cycle.

The devices are for use in weighing precious metals and other commodities of comparable value, or determination of dockage in the grain industry, and must be installed on a solid, stable surface free from vibration and in an environment free from air currents.

The devices are exempt from the requirements of sections 8, 9.2 and 9.6 of SGM3.

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APPROVAL:


The design, composition, construction and performance of the device type(s) identified herein have been evaluated in accordance with regulations and specifications established under the Weights and Measures Act. Approval is granted accordingly, pursuant to subsection 3(1) of the said Act.

The marking, installation, use and manner of use in trade of devices are subject to inspection in accordance with regulations and specifications established under the Weights and Measures Act. Certification of conformity is required in addition to this approval. Requirements relating to marking are set forth in sections 18 to 26 of the Weights and Measures Regulations. Requirements relating to installation, use and manner of use are set forth in Part V and in specifications established pursuant to section 27 of the said Regulations. Inquiries regarding inspection and certification of conformity should be addressed to the local inspection office of Consumer and Corporate Affairs Canada.

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 W.R. Virtue

Chief
Legal Metrology Laboratories

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