MAY 9 1990  
M A I

**NOTICE OF APPROVAL**

Issued by statutory authority of the Minister of Consumer and Corporate Affairs Canada for (category of device):

Electronic Top Loading Balance

**APPLICANT / REQUÉRANT:**

Fisher Scientific Limited  
112 Colonnade Road  
Nepean, Ontario  
K2E 7L6

**MODEL(S) / MODÈLE(S):**

CM150  
CM200  
PM100  
PM200  
PM300  
PM400  
PM460  
PM480  
PM600  
PM2000  
PM3000  
PM4000  
PM4600  
PM4800  
PM6000  
PM6

**NOTE:** This approval applies only to devices, the design, composition, construction and performance of which are, in every material respect, identical to that described in the material submitted, and that are typified by samples submitted by the applicant for evaluation for approval in accordance with sections 14 and 15 of the Weights and Measures Regulations. The following is a summary of principal features only.

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**SUMMARY DESCRIPTION:**

The approved devices are top loading electronic balances. Weight is sensed using an electromagnetic force compensation principle. All models are fitted with lockable levelling feet and a bullseye level.

The devices are equipped with a single operator control which performs the following functions:

To switch ON - press bar.

To switch OFF - lift bar.

To Zero - press bar. It operates to full scale capacity.

To Tare - press bar. It operates to full scale capacity.

A graphic display "Delta Trac" on the right side of the weight display gives an idea of the used and unused portions of the weighing range. When the scale is used, the Delta Trac becomes a plus(+) or minus (-) indicator showing how much weight is on the scale in relation to a target weight.

Blanking for over capacity is indicated by a series of flat bars on the weight display.

Delta Range - this is a mode that decreases the increment size by a factor of 10 and reduces capacity.

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SUMMARY DESCRIPTION: Continued

DESCRIPTION SOMMAIRE: Suite

The models PM and CM incorporate a Delta Display whereby the last decimal place is blanked out on initial application of weight but as the weight settles it re-appears. This function speeds up the weighing cycle.

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The devices are for use in weighing precious metals and other commodities of comparable value, or determination of dockage in the grain industry, or for determining the net content of packages that are not weighed in the presence of the purchaser at the time of sale.

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These devices must be installed on a solid, stable surface free from vibration and in an environment free from air currents.

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The devices are exempt from the requirements of sections 8, 9.2 and 9.6 of SGM3.

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NOTE:

On late models the Company name and model N° may be on the operator control bar.

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
## APPROVAL:

The design, composition, construction and performance of the device type(s) identified herein have been evaluated in accordance with regulations and specifications established under the Weights and Measures Act. Approval is hereby granted accordingly pursuant to subsection 3(1) of the said Act.

The marking, installation, use and manner of use in trade of devices are subject to inspection in accordance with regulations and specifications established under the Weights and Measures Act. Requirements relating to marking are set forth in sections 18 to 26 of the Weights and Measures Regulations. Installation and use requirements are set forth in Part V and in specifications established pursuant to section 27 of the said Regulations. A verification of conformity is required in addition to this approval. Inquiries regarding inspection and verification should be addressed to the local inspection office of Consumer and Corporate Affairs Canada.

W.R. Virtue

Chief,  
Legal Metrology Laboratories

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