NOV 17 1994

**NOTICE OF APPROVAL**

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Issued by statutory authority of the Minister of Industry,  
Science and Technology for:

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**NOTE:** This approval applies only to devices, the design, composition, construction and performance of which are, in every material respect, identical to that described in the material submitted, and that are typified by samples submitted by the applicant for evaluation for approval in accordance with sections 14 and 15 of the Weights and Measures Regulations. The following is a summary of the principal features only.

### SUMMARY DESCRIPTION:

The approved devices are top loading electronic balances. Weight is sensed using an electromagnetic force compensation principle. All models are fitted with lockable levelling feet and a bullseye level.

The devices are equipped with a single operator control which performs the following functions:

To switch ON - press bar.

To switch OFF - lift bar.

To Zero - press bar. It operates to full scale capacity.

To Tare - press bar. It operates to full scale capacity.

A graphic display "Delta Trac" on the right side of the weight display gives an idea of the used and unused portions of the weighing range. When the scale is used, the Delta Trac becomes a plus(+) or minus (-) indicator showing how much weight is on the scale in relation to a target weight.

Blanking for over capacity is indicated by a series of flat bars on the weight display.

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**SUMMARY DESCRIPTION: Cont'd**

Delta Range - this is a mode that decreases the increment size by a factor of 10 and reduces capacity.

The models PM and CM incorporate a Delta Display whereby the last decimal place is blanked out on initial application of weight but as the weight settles it re-appears. This function speeds up the weighing cycle.

The devices are for use in weighing precious metals and other commodities of comparable value, or determination of dockage in the grain industry, or for determining the net content of packages that are not weighed in the presence of the purchaser at the time of sale.

These devices must be installed on a solid, stable surface free from vibration and in an environment free from air currents.

The devices are exempt from the requirements of sections 8, 9.2 and 9.6 of SGM3.

<u>Model Number</u> <u>N° de modèle</u>	<u>Capacity</u> <u>Capacité</u>	<u>Increment</u> <u>Graduation</u>	<u>Delta (x10) Range</u> <u>Échelle Delta (x10)</u>
CM150	110g/150 carat	0.01g/0.001 carat	
CM200	110g/200 carat	0.01g/0.001 carat	
PM100	110g/550 carat	0.001g/0.01 carat/0.00001 lb	
PM200	210g/1050.45 carat	1mg/0.01 carat/0.00001 lb	
PM400	410g/2050.45 carat	1mg/0.01 carat/0.00001 lb	
PM460	410g/2050.4 carat	10mg/0.1 carat	60g/300 carat 1mg/0.01 ct.
PM480	410g/2050.4 carat	10mg/0.1 carat/0.0001 lb	80g/400 carat 1mg/0.01 ct. OR/ou 0.177 lb/ 0.00001 lb
PM600	610g/3054.5 carat	10mg/0.1 carat	
PM2000	2100g/10504.5 carat	10mg/0.1 carat/0.001 lb	
PM3000	3100g/15504 carat	100mg/1 carat/0.001 lb	
PM4000	4100g/20504.5 carat	10mg/0.1 carat/0.0001 lb	
PM4600	4100g/20504 carat	100mg/1 carat	600g/3000 ct. 10mg/0.1 ct.
PM4800	4100g/20504 carat	100mg/1 carat/0.001 lb	800g/4000 ct. 10mg/0.1 ct. OR/ou 1.764 lb/0.0001 lb
PM6000	6100g/20504 carat	100mg/0.001 lb	
PM6	6.1 kg	1g/0.01 lb	

**NOTE:**

On more recent models the Company name and model number may be on the operator control bar.

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**SUMMARY DESCRIPTION: Cont'd**

Revision 4 is to add the avoirdupois unit of measure and to add that the balance functions may be active.

**EVALUATED BY:**



Revision 4, Carl Cotton  
Calibration and Approvals Technologist  
Tel: (613) 952-0641

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

CM150 and/et CM200

**APPROVAL:**

The design, composition, construction and performance of the device type(s) identified herein have been evaluated in accordance with regulations and specifications established under the Weights and Measures Act. Approval is hereby granted accordingly pursuant to subsection 3(1) of the said Act.

The marking, installation and manner of use of trade devices are subject to inspection in accordance with regulations and specifications established under the Weights and Measures Act. Requirements relating to marking are set forth in sections 18 to 26 of the Weights and Measures Regulations. Installation and use requirements are set forth in Part V and in specifications established pursuant to section 27 of the said Regulations. A verification of conformity is required in addition to this approval. Inquiries regarding inspection and verification should be addressed to the local inspection office of Industry Canada.

D. W. Morgan

Manager,  
Weights and Measures Laboratories

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