Consumer and
Corporate Affairs Canada

Legal Metrology

Consommation
et Corporations Canada

Métrieologie légele

APPROVAL No. - N° D'APPROBATION

S.WA-4095 REV. 1

JAN 24 1989

NOTICE OF APPROVAL

Issued by statutory authority of the Minister of Consumer and Corporate Affairs Canada for (category of device):

Electro-Mechanical Precision Balance

APPLICANT / REQUÉRANT:

Sartorius GmbH
Weender Landstrasse 94-108
D-3400 Gottingen
Federal Republic of Germany

MODEL(S) / MODÈLE(S):

G150D-*K2

NOTE: This approval applies only to devices, the design, composition, construction and performance of which are, in every material respect, identical to that described in the material submitted, and that are typified by samples submitted by the applicant for evaluation for approval in accordance with sections 14 and 15 of the Weights and Measures Regulations. The following is a summary of principal features only.

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SUMMARY DESCRIPTION:

The approved device is a top loading electro-mechanical balance. Weight is sensed using an electro-magnetic force compensation principle. The balance is fitted with lockable levelling feet and a bullseye level and must be operated with a draft shield as the weight increment and sensitivity are such that air currents will prevent the balance from having a true zero. This balance is fitted with a + and - indication to the left of the weight display; when both are off, a zero display indicates that the zero position is accurate to $\pm 1/4$ of an increment.

The operating programs for this device are coded in a non-volatile memory (EEPROM Electronically Erasable PROM). These codes are accessible for change through the keyboard once an internal switch has been enabled. When the device has dual units of measure capability, there will be two codes for calibration, both will be disabled to meet the requirements of SGM3.10. This will also disable the "CAL" switch.

The device is equipped with operator controls which perform the following functions:

ON/OFF - Switches the display ON or OFF.

CAL - to calibrate the balance. An internal switch locks this system.

W1/W2 - Switches the display from carats to grams or vice versa.

TARE - Returns weight of empty container to zero balance. This key may also be used to switch display ON.

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Three holes have been machined in the base. Two in the front face and one underneath. These holes provide access to the calibration lock, linearity adjustment, and temperature resistor adjustment.

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
This device is for use in weighing precious metals and other commodities of comparable value, or determination of dockage in the grain industry, and must be installed on a solid, stable surface free from vibration and in an environment free from air currents.

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This device is exempt from the requirements of sections 8, 9.2 and 9.6 of SGM3.

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APPROVAL:

The design, composition, construction and performance of the device type(s) identified herein have been evaluated in accordance with regulations and specifications established under the Weights and Measures Act. Approval is hereby granted accordingly pursuant to subsection 3(1) of the said Act.

The marking, installation, use and manner of use in trade of devices are subject to inspection in accordance with regulations and specifications established under the Weights and Measures Act. Requirements relating to marking are set forth in sections 18 to 26 of the Weights and Measures Regulations. Installation and use requirements are set forth in Part V and in specifications established pursuant to section 27 of the said Regulations. A verification of conformity is required in addition to this approval. Inquiries regarding inspection and verification should be addressed to the local inspection office of Consumer and Corporate Affairs Canada.

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W.R. Virtue

Chief,
Legal Metrology Laboratories

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