**NOTICE OF APPROVAL**

Issued by statutory authority of the Minister of  
Industry for:

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**RATING / CLASSEMENT**

15°C to/à 35 °C

**NOTE:** This approval applies only to meters, the design, composition, construction and performance of which are, in every material respect, identical to that described in the material submitted, and that are typified by samples submitted by the applicant for evaluation for approval in accordance with sections 13 and 14 of the *Electricity and Gas Inspection Regulations*. The following is a summary of the principal features only.

### SUMMARY DESCRIPTION :

The Agilent 6890N (G1540N) Series Gas Chromatograph analyzes the composition of natural gas. The gas chromatograph separates the following components of natural gas:

- hydrogen
- helium
- oxygen
- nitrogen
- methane
- carbon dioxide
- ethane
- propane
- iso-butane
- n-butane
- iso-pentane
- n-pentane
- hexanes
- heptanes
- octanes
- nonanes
- decanes

The individual component concentrations are used to calculate the relative density and heating value of the sample. The detection of components is performed by a TCD (thermal conductivity detector) and an FID (flame ionization detector) which provide an output signal proportional to the concentration.

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Natural gas is injected into the sampling valve and is carried by the carrier gasses, helium and argon through five columns, inside the oven. The different components move through the columns at different rates, thereby allowing identification based on the retention time of each component.

The operator interface consists of a display, status board, and keyboard on the face of the gas chromatograph as well as a user interface program (EZ Chrom software) running on a personal computer. The EZ Chrom software analyses the data and displays the chromatogram. The EZ Chrom software version number can be viewed by selecting the "Help" menu and then "About" in the EZ Chrom software program.

Two Atco software applications called DGIS (Distribution Gas Information System) and Rumba, GMS (Gas Management System) are used to perform calibrations and analysis using GPA2261-00 and ASTM D1945 and to determine heating value and relative density at 15°C and 101.325 kPa using the GPA 2172 procedure and the GPA 2145-03 factors. The program version can be viewed when the program is opened.

The Agilent 6890N gas chromatograph's firmware version can be viewed by powering down the gas chromatograph and then powering up again. The firmware version will appear on the screen during boot up.

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or

10' x 1/8" - SS Mole Sieve  
13x45/80

Column 4: 4' x 1/8" - SS Haysep Q  
80/100, and  
4' x 1/8" - SS Mole Sieve  
13x45/60

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### Marking Requirements

Marking shall be in accordance with section 3.5 of LMB-EG-08 and section 20-3.1 of Provisional Specification S-G-03.

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### Modification Acceptance Letters (MALs)

MAL-G257

The purpose is to add new software and firmware versions and make corrections to some of the versions identified in the NOA. In addition to firmware version N.05.04 identified in the Notice of Approval, version N.05.05 is also approved for use.

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The software version 3.1.4 for the EZ Chrom Software identified in the Notice of Approval is obsolete and has never been put into production. The software version is now 3.1.7. Similarly, the software version 2.0 in the approval for the DGIS is also obsolete and is now 2.0.3056.29181. The Atco Rumba, GMS software identified in the approval has never been employed and is obsolete.

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#### MAL-G257 Rev.1-8

The purpose is to add the following DGIS software versions :

- 2009-05-26 – Version 2.0.3396.24262
- 2011-01-24 – Version 2.0.4038.12075
- 2011-06-16 – Version 2.0.4134.12164
- 2011-07-25 – Version 2.0.4219.10631
- 2012-05-02 – Version 2.0.4505.21683
- 2013-05-23 - Version 2.0.4750.20421
- 2013-07-10 – Version 2.0.4917.27186

In addition, on 2009-07-06 software version 3.3.2 was added for the EZ Chrom Software.

#### Revisions

Revision 1 added DGIS software version 2.0.5374.1180, integrated previous MALs, and performed minor editorial changes.

#### Evaluated By

AG-0529  
Judy Farwick  
Senior Legal Metrologist

AG-0529, Rev.1  
Allan Faust  
Junior Legal Metrologist

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**AG-0529 Rev.1**

## APPROVAL:

The design, composition, construction and performance of the meter type(s) identified herein have been evaluated in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. Approval is hereby granted accordingly pursuant to subsection 9(4) of the said Act.

The sealing, marking, installation, use and manner of use of meters are subject to inspection in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. The sealing and marking requirements are set forth in specifications established pursuant to section 18 of the *Electricity and Gas Inspection Regulations*. Installation and use requirements are set forth in specifications established pursuant to section 12 of the *Regulations*. Verification of conformity is required in addition to this approval for all metering devices excepting instrument transformers. Inquiries regarding inspection and verification should be addressed to the local office of Measurement Canada.

Patrick J. Hardock, P.Eng.  
Senior Engineer – Gas Measurement  
Engineering and Laboratory Services Directorate

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AG-0529, Rév. 1