NOTICE OF APPROVAL

Issued by statutory authority of the Minister of
Industry for:

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RATING / CLASSEMENT

15°C to/à 35 °C

NOTE: This approval applies only to meters, the design, composition, construction and performance of which are, in every material respect, identical to that described in the material submitted, and that are typified by samples submitted by the applicant for evaluation for approval in accordance with sections 13 and 14 of the *Electricity and Gas Inspection Regulations*. The following is a summary of the principal features only.

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SUMMARY DESCRIPTION :

DESCRIPTION SOMMAIRE

The Agilent 6890N (G1540N) Series Gas Chromatograph analyzes the composition of natural gas. The gas chromatograph separates the following components of natural gas:

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The individual component concentrations are used to calculate the relative density and heating value of the sample. The detection of components is performed by a TCD (thermal conductivity detector) and an FID (flame ionization detector) which provide an output signal proportional to the concentration.

Interfacing with the gas chromatograph is performed by using the keyboard and display located on the front of the device. For OpenLab CDS 2.4 the interfacing is done using an AIC – Agilent Instrument Controller.

The gas chromatograph is also connected to a computer which uses EZ Chrom or OpenLab CDS 2.4 software to analyse the gas composition. Calibration of the gas chromatograph is performed through the EZ Chrom or OpenLab CDS 2.4 software's using ASTM D1945-14 and GPA2261-13 procedure.

ATCO software application, DGIS (Distribution Gas Information System) is used in the analysis to determine real gross heating value and relative density at 15°C and 101.325 kPa (dry air, dry gas). The GPA 2172-14 calculations and the GPA 2145-16 factors are used in determining the heating value and relative density.

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Column 1: 3' x 1/8" - SS Haysep Q
- 80/100

Column 2: 4 x 1/8" - SS Haysep N
- 80/100
or
6' x 1/8" Haysep N 80/100

Column 3: 6' x 1/8" - SS Mole Sieve
13x45/60
or
10' x 1/8" - SS Mole Sieve
13x45/60
or
10' x 1/8" - SS Mole Sieve
13x45/80

Column 4: 4' x 1/8" - SS Haysep Q
80/100, and
4' x 1/8" - SS Mole Sieve
13x45/60
or
3' x 1/8" Haysep Q 80/100,
and 4' x 1/8" Molsieve 5A 60/80

Column 5: HP1 50m x 0.32 mm x
1.05 :m film

Carrier Gas: Helium, Argon and Hydrogen

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Microprogramme: N.05.04, N.05.05, N.05.06,
N.06.07

Marking Requirements

The Agilent 6890N gas chromatograph's firmware version can be viewed by navigating to *Options* → *Diagnostic* → *Instrument Status* → *Version* on the front interface of the gas chromatograph

The EZ Chrom software version number can be viewed by selecting the "Help" menu and then "About" in the EZ Chrom software program.

The OpenLab CDS 2.4 software version number can be viewed in the configuration settings of the control panel.

The DGIS program version can be viewed when the program is opened.

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Photographs and Diagrams / Photos et diagrammes


Figure 1: Typical Nameplate / Plaque signalétique typique

Figure 2: Agilent 6890N

APPROVAL:

The design, composition, construction and performance of the meter type(s) identified herein have been evaluated in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. Approval is hereby granted accordingly pursuant to subsection 9(4) of the said Act.

The sealing, marking, installation, use and manner of use of meters are subject to inspection in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. The sealing and marking requirements are set forth in specifications established pursuant to section 18 of the *Electricity and Gas Inspection Regulations*. Installation and use requirements are set forth in specifications established pursuant to section 12 of the *Regulations*. Verification of conformity is required in addition to this approval for all metering devices excepting instrument transformers. Inquiries regarding inspection and verification should be addressed to the local office of Measurement Canada.

Original copy signed by:

Christian Lachance, P.Eng.
Senior Engineer – Gas Measurement
Engineering and Laboratory Services Directorate

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