**NOTICE OF APPROVAL**

Issued by statutory authority of the Minister of  
Industry for:

**TYPE OF DEVICE**

Gas Chromatograph

**APPLICANT**

Scion Instruments Canada Limited  
Unit 2008 – 63rd Avenue East  
Edmonton, Alberta  
T6E 0G9

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**NOTE:** This approval applies only to meters, the design, composition, construction and performance of which are, in every material respect, identical to that described in the material submitted, and that are typified by samples submitted by the applicant for evaluation for approval in accordance with sections 13 and 14 of the *Electricity and Gas Inspection Regulations*. The following is a summary of the principal features only.

### SUMMARY DESCRIPTION :

The Bruker 450 and Scion 456GC gas chromatographs use two TCD and one FID detectors. Helium and Argon are used as carrier gases for the TCD detector and Hydrogen gas is used with the FID detector. Individual component concentration is determined by the chromatograph with the use of a single injection of the sample gas. This sample fills three sample loops. The front TCD detector uses a 12 foot mole sieve column with an Argon carrier gas to separate He, H<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub> and C<sub>1</sub>. The middle TCD detector uses a 6 foot Haysep N 80/100 mesh column with a Helium carrier gas to separate CO<sub>2</sub> and C<sub>2</sub>. The FID detector uses 10m X 0.15 mm CPSIL 5 CB capillary column with a Hydrogen carrier gas to separate C<sub>3</sub> to C<sub>10+</sub>.

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- 120 VAC, 50/60 Hz, 2400 VA

### Carrier Gas:

- Helium for column #1
- Argon for column #2
- Hydrogen for column #3

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### Column Information:

- Column #1: Hayesep R (2' x 1/8" Ultimetmetal 80/100 Mesh) connected to and before a Molesieve 5 A (4 m x 1/8" Ultimetmetal 80/100 Mesh) connected to a TCD.
- Column #2: Hayesep N (6' x 1/8" Ultimetmetal 80/100 Mesh) connected to a TCD.
- Column #3: CP-Sil 5CB (10 m x 0.15 mm with 2 µm film thickness) connected to a FID.

### Gas Composition

The Bruker 450 and Scion 456GC are approved for determining relative concentration of the following components in a gas sample:

- N<sub>2</sub>, nitrogen
- O<sub>2</sub>, oxygen
- CO<sub>2</sub>, carbon dioxide
- CH<sub>4</sub>, methane
- C<sub>2</sub>H<sub>6</sub>, ethane
- C<sub>3</sub>H<sub>8</sub>, propane
- i-C<sub>4</sub>H<sub>10</sub>, iso-butane
- n-C<sub>4</sub>H<sub>10</sub>, normal butane
- i-C<sub>5</sub>H<sub>12</sub>, iso-pentane
- n-C<sub>5</sub>H<sub>12</sub>, normal pentane
- C<sub>6</sub>H<sub>14+</sub>, hexane and above

### Energy Calculations:

The SGS Canada Inc. Excel spreadsheet uses the outlines in GPA 2286 for determining the heating value and relative density from the sample gas composition determined by the chromatograph.

### Software:

- Instrument Firmware: Versions 3.03, 4, and 5.08
- Workstation Software: Galaxie Version 1.9, and Compass Versions 3.0 and 3.01

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## MARKING REQUIREMENTS

Markings shall be in accordance with Section 3.5 and 20-3.1 of LMB-EG-08 and section 20-3.1 of S-G-03 (2008-10-08).

Displayed is on the printed report:  
Base temperature and pressure for the calorific power.

## VERIFICATION

Energy Calculation

The true calorific power and relative density of the sample gas shall be calculated using the true relative concentrations of all components of the sample gas as established in accordance with the Gas Processors Association standards, *Calculation of Gross Heating Value, Relative Density and Compressibility Factor for Natural Gas Mixture from Compositional Analysis*, GPA 2172, and *Table of Physical Constants for Hydrocarbons and Other Compounds of Interest to the Natural Industry*, GPA 2145.

## REVISIONS

MAL-G319:

MAL-G319 approves the instrument firmware version 4, and workstation software Compass Version 3.0.

Revision 1:

Revision 1 approves the Scion 456GC (re-branding) gas chromatograph, and instrument firmware version 5.08.

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## EVALUATED BY

### Original:

Claude Dupont  
Senior Legal Metrologist

### MAL-G319:

Randy Byrtus  
Manager, Gas Measurement

### Revision 1:

Adam Falicki  
Junior Legal Metrologist

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**AG-0591 Rev. 1**

## APPROVAL:

The design, composition, construction and performance of the meter type(s) identified herein have been evaluated in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. Approval is hereby granted accordingly pursuant to subsection 9(4) of the said Act.

The sealing, marking, installation, use and manner of use of meters are subject to inspection in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. The sealing and marking requirements are set forth in specifications established pursuant to section 18 of the *Electricity and Gas Inspection Regulations*. Installation and use requirements are set forth in specifications established pursuant to section 12 of the Regulations. Verification of conformity is required in addition to this approval for all metering devices excepting instrument transformers. Inquiries regarding inspection and verification should be addressed to the local office of Measurement Canada.

## Original copy signed by:

Christian Lachance, P.Eng.  
Senior Engineer – Gas Measurement  
Engineering and Laboratory Services Directorate

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