NOTICE OF APPROVAL

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MANUFACTURER

FABRICANT

Bruker Chemical Analysis Division
555 Steeles Avenue East
Milton, Ontario, Canada
L9T 1Y6

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NOTE: This approval applies only to meters, the design, composition, construction and performance of which are, in every material respect, identical to that described in the material submitted, and that are typified by samples submitted by the applicant for evaluation for approval in accordance with sections 13 and 14 of the *Electricity and Gas Inspection Regulations*. The following is a summary of the principal features only.

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SUMMARY DESCRIPTION :

DESCRIPTION SOMMAIRE :

The Varian 450-GC gas chromatograph uses two TCD and one FID detectors. Helium, Hydrogen or Argon is used as carrier gases for the TCD detector and Hydrogen gas is used with the FID detector. Individual component concentration is determined by the chromatograph with the use of a single injection of the sample gas. This sample fills three sample loops. The front TCD detector uses a nickel 12 foot x 1/8" mole sieve column with an Argon carrier gas to separate He, H₂, O₂, N₂ and C₁. The middle TCD detector uses a nickel 6 foot x 1/8" Haysep N 80/100 mesh column with a Helium or Hydrogen carrier gas to separate CO₂ and C₂. The FID detector uses 10m X 0.15mm CP-SIL 5 CB capillary column with a Hydrogen carrier gas to separate C₃ to C₁₅₊.

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Column Information:

- Column #1: Haysep N (6' x 1/8" Nickel 80/100 mesh)
- Column #2: Mole Sieve 13x packed, nickel 12' x 1/8" 45/60 mesh.
- Column #3: CP-Sil 5CB (10m x 0.15mm with 2µm film thickness) connected to a FID.

Gas Composition

The 450-GC identifies gas components up to C₁₅₊ and sums C₁₀₊ components for the purpose of calculating heating value and density.

Energy Calculations:

The software program Norlims Sample Processing is used for determining the air correction, heating value and relative density from the sample gas composition determined by the chromatograph in accordance with GPA2172 and AGA 5.

Software / Firmware

- Instrument Firmware: version 3.02, 4.05
- Workstation Software: Galaxie version 1.9.302.952 or Compass CDS: version 3.0.0.68
- Norlims Sample Processing: version 5.12, 5.17

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Marking Requirements

Markings shall be in accordance with Section 3.5 and 20-3.1 of LMB-EG-08 and section 20-3.1 of S-G-03 (2008-10-08).

Displayed is on the printed report: Base temperature and pressure for the calorific power.

VERIFICATION

The true calorific power and relative density of the sample gas shall be calculated using the true relative concentrations of all components of the sample gas as established in accordance with the Gas Processors Association standards, *Calculation of Gross Heating Value, Relative Density and Compressibility Factor for Natural Gas Mixture from Compositional Analysis*, GPA 2172, and *Table of Physical Constants for Hydrocarbons and Other Compounds of Interest to the Natural Industry*, GPA 2145.

REVISION AND MAL'S

Revision 1:
Revision 1 add the Compass CDS analytical software application version 3.0.0.68, add a version 5.17 of the NorLIMS analytical software, add a new firmware version 4.05 and allow for either helium or hydrogen to be the carrier gas for column #1, middle TCD.

EVALUATED BY

Original: (2011-11-02)
Randy Byrtus
Manager, Gas Measurement

Revision 1:
Claude Dupont
Senior Legal Metrologist

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Figure 1- Chromatograph / Chromatographe

Figure 2 - column / colonne

Measurement
Canada

An Agency of
Industry Canada

Mesures
Canada

Un organisme
d'Industrie Canada

APPROVAL No. - N° D'APPROBATION

AG-0598 Rev 1

APPROVAL:

The design, composition, construction and performance of the meter type(s) identified herein have been evaluated in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. Approval is hereby granted accordingly pursuant to subsection 9(4) of the said Act.

The sealing, marking, installation, use and manner of use of meters are subject to inspection in accordance with regulations and specifications established under the *Electricity and Gas Inspection Act*. The sealing and marking requirements are set forth in specifications established pursuant to section 18 of the Electricity and Gas Inspection Regulations. Installation and use requirements are set forth in specifications established pursuant to section 12 of the Regulations. Verification of conformity is required in addition to this approval for all metering devices excepting instrument transformers. Inquiries regarding inspection and verification should be addressed to the local office of Measurement Canada.

Original copy signed by :

Christian Lachance, P.Eng.
Senior Engineer – Gas Measurement
Engineering and Laboratory Services Directorate

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